

# Staebel P5500 Homecare Bed

Product Code: 7-P5500

The Staebel P5500 Homecare Bed is designed for easy transportation and storage.





The Staebel P5500 Homecare Bed is designed for easy transportation and storage. Ideal for short term cases, the P5500 bed can be assembled and dismantled with ease through the incorporation of a simple lock system which is secured with a pin, and folds compactly for storage.

### Product Features:

- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 200kg
- Back Lit Handset
- Ironstone Powder Coated Frame



Forward & Reverse Trendelenburg



Chair Position Function



100mm Braking Castors



Bluetooth Connectivity  
\*\*\*handset must be ordered separately\*\*\*



### Accessories – Added after manufacture

- Mattresses
- Bed Extension
- Horizontal Bed Rails
- Clamp on Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Pole
- UPS Emergency Power Bed Unit
- IV Pole
- Combo Grap Handle
- Reading Lamp
- Battery Backup



### Specifications:

Size	Single	King Single
External Dimensions	930 x 2190 mm	1090 x 2190 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm
Hi-Lo Position	295 – 700 mm	295 – 700 mm
Back / Knee Angle	66° / 35°	66° / 35°
TR / ATR Angle	11° / 11°	11° / 11°
SWL	200 kg	200 kg