

Alrick Plexa Bed

Product Code: 7-PLEXA

Designed to provide safe and comfortable sleeping arrangements for users with complex needs. This bed includes full length drop down rails to provide security. It features a locking system, ensuring the rails stay in place while the user is sleeping.



Product Features:

- Australian Made and Owned
- Floorline Model
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 250kg SWL
- Full Length Drop Down Rails
- Back Lit Handset
- Square Natural Oak End Boards
- Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



Chair Position
Function



75mm Twin
Wheel Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Specifications:

Size	Single	King Single	Double
External Dimensions	1040 x 2150 mm	1210 x 2150 mm	1510 x 2150 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm	1370 x 2000 mm
Hi-Lo Position	150-725 mm	150-725 mm	170-725 mm
Rail Hi-Lo Height	100-400 mm	100-400 mm	100-400 mm
Backrest / Knee Angle	62° / 35°	62° / 35°	62° / 35°
TR / ATR Angle	15° / 15°	15° / 15°	15° / 15°
SWL	250 kg	250 kg	250 kg

Issue Date: 130224
Revision Date: 240426 | V5



Providing a safe & comfortable sleeping arrangement for users with complex needs.

Optional Features – Added at time of manufacture

- Dual Locking Castors
- Bluetooth Handset (HC40)
- Battery Back Up
- Under Bed Light (Warm White)
- Under Bed Light with Custom Colours
- USB Charging Ports
- Custom End Board Styles and Colours

Custom options are available on request at an additional cost (\$)
– Contact the team for further details

Stocked Endboard Colours & Styles

Custom Endboard, styles and colours also available.



Laminate
Square Headboard & Footboard




Upholstered
Plain Headboard & Footboard



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered

Accessories – Added after manufacture

- Mattresses
- Self Help Pole
- Emergency Bed Power Unit
- Bluetooth Handset (7-BT-4-HB200)
- IV Pole
- EN T-Bar Buffer
- Reading Lamp
- Padded Side Protectors



Custom Jersey
Endboards



Padded Side
Protectors

