

Alrick Zyra Homecare Bed

Product Code: 7-ZYRA

The Alrick Zyra Homecare Bed is engineered with a robust design to ensure safety and stability. The bed's sturdy construction, reinforced frame, and high-quality materials contribute to its exceptional weight capacity, boasting an impressive safe working load capacity of 500kg. This bed delivers strength and durability from King Single width to King Double.




Product Features:

- Australian Made and Owned
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 500kg SWL
- Under Bed Light (Warm White)
- Pull Out Extension
- Back Lit Handset
- Battery Back Up
- Curved Natural Oak End Boards
- Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



Chair
Position



125mm Braking
& Directional
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***

Optional Features - Added at time of manufacture

- CPR Backrest
- Electric Vascular Support
- Central Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- USB Charging Ports
- Under Bed Light with Custom Colours
- Custom Powder Coat Colour
- Custom End Board Style and Colour
- Custom Length Sleeping Surface (2m+)

Accessories - Added after manufacture

- Mattresses
- Horizontal 400HR Bed Rails (pair)
- Padded Cotside Protectors (pair)
- Clamp On Support Rail
- StandEzy Straight Handle
- StandEzy Hook Handle
- Self Help Pole
- Bluetooth Handset (7-BT-5-HB200)
- 200mm Foam Bolster
- UPS Emergency Power Bed Unit
- IV Pole
- Combo Grab Handle
- T-Bar Buffer
- Reading Lamp

Ensuring safety and stability in a home care environment.

Stocked Endboard Styles & Colours

Custom Endboard styles & colours available upon request.



Laminate
Curved Style



Upholstered
Jersey Headboard &
Plain Footboard




Timber
Eltham Style




Natural Oak
Laminate



Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber



Saddle
Timber



Specifications:

Size	King Single	Extra Wide	Double	Queen Double	King Double
External Dimensions	1090 x 2050 mm	1210 x 2050 mm	1390 x 2050 mm	1520 x 2050 mm	1820 x 2050 mm
Sleeping Surface	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm	1800 x 2000 mm
Hi-Lo Position	360 - 800 mm	360 - 800 mm	360 - 800 mm	360 - 800 mm	360 - 800 mm
Underbed Clearance	160 mm	160 mm	160 mm	160 mm	160 mm
Back / Knee Angle	65° / 40°	65° / 40°	65° / 40°	65° / 40°	65° / 40°
TR / ATR Angle	16° / 16°	16° / 16°	16° / 16°	16° / 16°	16° / 16°
SWL	500 kg	500 kg	500 kg	500 kg	500 kg