

Customise your Alrick Floorline Bed

Prescription Form

Please fill out the below form to customise your Alrick Floorline Bed. Once completed, return the form to the team at Active Medical.



1. Form Information

Client Information

Full Name				Contact Name			
Street Address							
Suburb				Contact Number/Email			
State		Postcode		Plan Date - Start to Finish		to	
Client Weight	kg	Funding Number		Plan Manager			

NDIS	Homecare Package	Other: _____
------	------------------	--------------

Prescriber Information

2. Sleeping Surface Measurements

900mm (Single)	1070mm (King Single)	1200mm (Extra Wide)
1370mm (Double)	1500mm (Queen Double)	1800mm (King Double)

Length



2000mm (Standard) (go to section 5) Custom Length (complete section 3)

3. Customise the Sleeping Surface

The bed is split into 4 separate sections including:

1. Backrest
2. Fixed
3. Small Kneebreak
4. Large Kneebreak

Please see standard measurements to the right. Customise each sleeping surface sections to suit user needs using diagram & measurements below.



Custom Length:

Backrest	mm	Fixed	mm	Small Kneebreak	mm	Large Kneebreak	mm
Total Length							mm

4. Bed Requirements

SWL Required





Safe Working Load is **200kg**.

To increase Safe Working Load, see 300kg Handset Upgrade (Single & King Single Only).







Castors

75mm Twin Wheel Braking (Standard)	Dual Locking (\$)	100mm Twin Wheel Braking (\$)
------------------------------------	-------------------	-------------------------------

Controls

Backlit Handset (Standard) 	Attendant Control (\$) 	Embedded Foot Control (\$)  <small>(Only available in Laminate or Eltham Endboards)</small>	HC40 Bluetooth Handset (\$) <small>Note: The standard handset will still be supplied</small> 
---	---	---	--

Optional Features

Electric Vascular Support (\$) 	CPR Backrest (\$) 	Battery Backup (\$) 
USB Charger Port (\$) <small>(Specify Quantity & which side of the bed)</small> 	Underbed Lighting (White) (\$) <small>(Specify Quantity & which side under the bed)</small> 	Underbed Lighting with Mood Lighting (all colours included) (\$) <small>(Specify Quantity & which side under the bed)</small> 

Specify Here: _____

5. Bed Frame Design

Powdercoat Colours

Ironstone (Standard)	APO Grey (\$)	Rivergum Beige (\$)	Signal Red (\$)	Manor Red (\$)	Orange X15 (\$)
Lemon Yellow (\$)	Mint Green (\$)	Cabana Green (\$)	French Blue (\$)	Blaze Blue (\$)	Navy (\$)
Dark Violet (\$)	Pink Blush (\$)	Precious Silver Pearl (\$)	Pearl White (\$)	Black Satin (\$)	

Other — (please specify colour & range)



6. Endboard Design

First, select your style of endboard then your preferred colour/stain.

Laminate (move to Laminate Section)	Timber (move to Timber Section)	Upholstered (move to Upholstered Section)
-------------------------------------	---------------------------------	---

Laminate - Endboards (Available as Head or Footboards)









The standard height above mattress dimensions includes: **Headboard is 440mm** above mattress base and **Square Low Footboard is 190mm** above mattress base. Customise each End Board Height to suit clients needs using below table.

Standard	Slender Arch	Curve	Queen Ann	Square
				
Slender Arch Headboard	Headboard	Headboard	Headboard	Headboard
Square Low Footboard	Footboard	Footboard	Footboard	Footboard

Custom End Board Height (\$)


Headboard Requested Height	mm	Footboard Requested Height	mm
----------------------------	----	----------------------------	----

Laminate Colours *Laminate colours from Polytec Wood Matt Range

Natural Oak (Standard)	Jamaican Walnut (\$)	South West Jarrah (\$)	European Walnut (\$)
			
Shannon Oak (\$)	Soft Walnut (\$)	Artisan Oak (\$)	Classic White (\$)
			









Cant find what you are after? — (please specify style & colour or attach image)

Timber Endboards

San Remo (\$)		Kensington (\$)	
Headboard		Headboard	
Footboard		Footboard	
Eltham (\$)		Rylestone (\$)	
Headboard		Headboard	
Footboard		Footboard	
Luxa Timber Surround (\$) (Headboard to be ordered separately)		Other — (please specify colour & range)	
Timber Foot Panel			
Side Panels (Pair)			

Cant find what you are after? — (please specify style & colour or attach image)

Timber Stain Colours *Colours / Stains only available for timber products

Beech (\$)	Burnt Walnut (\$)	Cognac (\$)	Chestnut (\$)
			
Natural (\$)	Saddle (\$)	Cocoa (\$)	Cinnamon (\$)
			









Cant find what you are after? — (please specify style & colour or attach image)

Upholstered Endboards

Loft (\$)			Princeton (\$)			Plain Rectangle (\$)		
	Headboard			Headboard			Headboard	
	Footboard			Footboard				
Brooklyn (\$)			Cambridge (\$)			Jersey		
	Headboard			Headboard			Headboard	
Panache (\$)			Plain Rectangle (\$)			Luxa Padded Surround (\$) (Headboard to be ordered separately)		
	Headboard			Footboard			Padded Foot Panel	
							Side Panels (Pair)	









Cant find what you are after? — (please specify style & colour or attach image)

Fabric Colours *Colours only available for upholstered products

Ebony (\$)	Ash (\$)	Almond (\$)	Ivory (\$)
			
Ocean (\$)	Opal (\$)	Dusk (\$)	Eucalyptus (\$)
			
















Cant find what you are after? — (please specify style & range)

Vinyl Colours *Colours only available for upholstered products

Black (\$)	Pewter (\$)	Brazil (\$)	Snow (\$)
			
Periwinkle (\$)	Emerald (\$)	Lemon Twist (\$)	Flames (\$)
			

Cant find what you are after? — (please specify style & colour or attach image)

Optional Accessories (\$)

 <p>StandEzy Hook Handle Bed Pole (7-EN9-SE-HH)</p>	 <p>Combo Grab Rail (7-CGH)</p>	 <p>Floorline T-Bar Buffer (7-EN9-TBR)</p>	 <p>I.V Pole (7-AIVP)</p>	 <p>Reading Light (7-ABSL)</p>
 <p>Padded Cotside Protectors (7-COTPRO-G)</p>	 <p>Horizontal Rails (7002HR-EN7(width)-IS)</p>	 <p>Horizontal Rails (7002HR-EN7(width)-APO)</p>	 <p>200mm Bed Extension (7-C.EXTEN9000-(width)-IS)</p>	 <p>200mm Bed Extension (7-C.EXTEN9000-(width)-APO)</p>
 <p>Foam Bolster (7-BLR)</p>	 <p>Self Help Pole (7-ASHP/H)</p>	 <p>300kg Handset Upgrade (Single & King Single Only)</p>	 <p>Bluetooth Handset (7-BT-4-HB200)</p>	 <p>UPS Emergency Power Supply (7-UPS-EBPU-NH)</p>

6. Please Note

All colours are an indication only and may vary. Whilst we will endeavour to accommodate all your requests, there may be some manufacturing limitations which will not allow us to provide all requests. Please use the Comment/Request box below.

(\$) – Indicates a change in price associated with the selected function, style, colour, or accessory.

7. Comments / Requests