

Nestie Nova Maternity Bed

Product Code: 7-NOVA

Designed with a wide height range, and removable leg stirrups, ensures comfort for bonding and feeding while supporting both users and care staff during post-birth care.



Specifications:

Size	King Single	Double	Queen Double
Overall Dimensions	1165 x 2040 mm	1465 x 2040 mm	1595 x 2040 mm
Sleeping Surface	1070 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Range	390 - 830 mm	390 - 830 mm	390 - 830 mm
Underbed Clearance	160 mm	160 mm	160 mm
Back / Knee Angle	77° / 40°	77° / 40°	77° / 40°
TR / ATR Angle	16° / 16°	16° / 16°	16° / 16°
SWL	300 kg	300 kg	300 kg

Issue Date: 010222
Revision Date: 270426 | V4



The Nestie Nova Maternity Bed has a high SWL of 300kg with height adjustability that empowers nurses to provide exceptional care and minimising strain. The Nestie Nova has optional removable leg stirrups which provides a safe, comfortable space, enabling post birth suturing and medical check-ups with ease. The Nestie Nova also features an adjustable backrest, enhancing the comfort while feeding and bonding.

Product Features:


- Australian Made and Owned
- Electric Hi-Low
- Electric Backrest
- Electric Kneebreak
- Battery Backup
- Underbed Light (Warm White)
- Backlit Handset
- Safe Working Load 300kg
- Ironstone steel frame and Matt Classic White Infill Endboards
- APO Grey and Ironstone Powdercoated Frame
- Birthing Stirrup Bar (*Stirrups to be purchased separately*)



Forward & Reverse
Trendelenburg



Manual CPR
Release



150mm Twin
Wheel Central
Locking Castors



Bluetooth
Connectivity
***handset must be
ordered separately***



Optional Features - Added at time of manufacture

- Upgrade to 400kg SWL
- Electric Vascular Support
- Underbed Light with Custom Colours
- Attendant Control
- Bluetooth Handset
- Laminate or Upholstered Endboards

Accessories - Added after manufacture

- Stirrups
- Mattresses
- IV Pole
- Self Help Pole & Triangle
- Bariatric Self Help Pole & Triangle
- 400HR Horizontal Deluxe Rails

