

Secura Mental Health Floorline Bed

Product Code: 7-SEC-FL

The Secura Mental Health Floorline bed provides a unique solution that combines both safety and clinical performance, ensuring that mental health facilities do not have to compromise one for the other.



Product Features:


- Australian Made
- Electric hi/low
- Electric kneebend
- Electric auto regression backrest
- Anti tamper screws
- Encased wiring and control box
- 75mm twin wheel breaking castors
- Bolt on head board 25mm thick (Natural Oak)
- Bolt on foot board 25mm thick with embedded bed controls (Natural Oak Standard)
- Ironstone powdercoated frame

Optional Features - Added at time of manufacture

- Under bed light
- Dual breaking Castors
- Trendelenburg
- Battery backup
- CPR Release
- Bluetooth Handset
- Custom endboard colours

Accessories - Added after time of manufacture

- Mattresses



Specifications:

Size	Single	King Single	Double	Queen Double
External Dimensions	920 x 2115 mm	1090 x 2115 mm	1390 x 2115 mm	1520 x 2115 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm
Hi-Lo Position	100 - 725 mm	100 - 725 mm	100 - 725 mm	100 - 725 mm
Backrest / Knee Angle	62° / 35°	62° / 35°	62° / 35°	62° / 35°
SWL	200 kg	200 kg	250 kg	250 kg