

Staebel Forearm Walker

Product Code: 680SFW

The Staebel Forearm Walker allows the user to move safely and independently during activities, such as walking, standing up, as well as sitting back down. It also serves as an indispensable rehabilitation aid when recovering from injuries.





Product Features:

- Height adjustable via gas spring with one simple action
- Horizontally adjustable handles
- Increase user space with adjustable armrests:
- Swivel castors with diameter of 100mm with brake
- Powdercoated metal frame in Light Grey
- Upholstered armrests in colour Antracit
- Sealed arm pads for infection control
- S Shaped IV Pole with 2 hooks (removeable)
- Oxygen Bottle Carrier (removeable)
- Soft Hand Brakes for comfortable patient control
- Available in Standard or Bariatric size

Accessories:

- IV Pole Holder with IV Pole (6820)
- Oxygen Bottle Holder (6821)

Specifications:

Product Code	6800	6808
Size	Standard	Bariatric
Overall Width	720 mm	870 mm
Width Between Arms (Internal Dimensions)	370 - 490 mm	550 - 680 mm
Width Between Arms (External Dimensions)	640 - 770 mm	820 - 950 mm
Height	930 - 1220 mm	1060 - 1360 mm
Overall Length	750 mm	750 mm
SWL	170 kg	300 kg