

Staebel Sovereign Wheelchair

Product Code: 4300C

The most universal wheelchair –
Adjustable to suit any user.



The Sovereign Wheelchair is designed to adapt effortlessly to the unique needs of each user, offering an exceptional level of versatility and comfort. Crafted with insights from occupational therapists, this wheelchair is built to meet a wide range of requirements, making it the ultimate all-in-one solution for clients.

Say goodbye to complicated customisations and endless add-ons! The Sovereign Wheelchair provides unmatched adjustability without the need for extra parts, seamlessly accommodating various user preferences. This one-chair-fits-all design makes it a reliable, ready-to-go option that simplifies patient support and enhances user experience.



125-170kg
SWL



Excellent User
Comfort



Lightweight
Construction



Foldable
Design

Product Features:

- Adjustable & removable legrests
- Flip back & removable armrests
- Adjustable armrest height & length
- Seat depth, height & angle adjustment
- Backrest height, tension & angle adjustment
- Crash Tested
- Safe Working Load ranges from 125kg to 170kg



Specifications:

Product Code	4340C - 37	4340C - 40	4340C - 43	4340C - 46	4340C - 49	4352C - 52	4352C - 55	4352C - 58
Seat Width	370 mm	400 mm	430 mm	460 mm	490 mm	520 mm	550 mm	580 mm
Overall Width	392 mm	422 mm	452 mm	482 mm	512 mm	542 mm	572 mm	602 mm
Anterior Height	420 - 500 mm	420 - 500 mm	420 - 500 mm	420 - 500 mm	420 - 500 mm	420 - 500 mm	420 - 500 mm	420 - 500 mm
Posterior Height	410 - 510 mm	410 - 510 mm	410 - 510 mm	410 - 510 mm	410 - 510 mm	410 - 510 mm	410 - 510 mm	410 - 510 mm
Seat Depth	400 - 500 mm	400 - 500 mm	400 - 500 mm	400 - 500 mm	400 - 500 mm	450 - 500 mm	450 - 500 mm	450 - 500 mm
Backrest Height (2cm increments)	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm
Backrest Angle (5° increments)	85° - 105°	85° - 105°	85° - 105°	85° - 105°	85° - 105°	85° - 105°	85° - 105°	85° - 105°
Lower Leg Length (2cm increments)	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm
Armrest Height (2cm increments)	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm
Armrest Length (4cm increments)	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm
Footplate Width	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm
NET/Gross Weight (Approximate)	17 kg / 21 kg	17 kg / 21 kg	18 kg / 21 kg	18 kg / 21 kg	18 kg / 21 kg	21 kg / 24 kg	21 kg / 24 kg	21 kg / 25 kg
SWL	125 kg	125 kg	125 kg	125 kg	125 kg	170 kg	170 kg	170 kg