

Staebe! Sovereign Wheelchair

Order Form

The Sovereign Wheelchair is a basic lightweight wheelchair that can be adjusted to suit almost many users.



Form Information

Full Name

Contact

Date

Wheelchair Specifications

Wheelchair (Safe Working Load 170kg)

Seat Width

370 mm

400 mm

430 mm

460 mm

490 mm

Seat Depth

400 mm

425 mm

450 mm

475 mm

500 mm

Wheelchair (Safe Working Load 170kg)

Seat Width

520 mm

550 mm

580 mm

Seat Depth

450 mm

475 mm

500 mm

Anterior Seat Height

420 mm

445 mm

450 mm

470 mm

475 mm

500 mm

Posterior Seat Height

410 mm

430 mm

450 mm

470 mm

490 mm

510 mm

Rear Wheel Specifications

Rear Wheel Size (* 12" transit option only available on products up to 490mm seat width – Not available on 170kg model)

24" Solid tyre – as standard

12" Transit option*

Rear Wheel Position


Position 1

Position 2

Position 3

Position 4

Specify position of rear wheel as per diagram. (B shows Rear Wheel Horizontal position).



Standard Hollow push Rim

Drum Brake Option, includes brake, cable and wheels

One Arm Drive Kit with double push rim

Mounted Left

Mounted Right

Rear Wheel Accessories

Amputee brackets, Pair

Anti tip bars with inbuilt tip assist, Pair

Brakes (Wheel Lock - Brake Lever Extension)

Left

Right

Backrest Height

400 mm

420 mm

440 mm

460 mm

Backrest Angle (If no specified, backrest angle standard at 90°)

85°

90°

95°

100°

105°

Footrest Setting- Footplate/lower leg length (Adjustable from 40-48cm - 2cm increments)

400 mm

420 mm

440 mm

460 mm

480 mm

Footplate Assembly Individual angle - adjustable footplates, adjustable to close gap)

Footplate width - adjustable 105-155mm. If applicable specify setting:

mm

Back Upholstery

Adaptable - upholstery colour black

Push Handles

Standard height adjustable and detachable

Side Panel

Flip Back Height - adjustable armrest with removable mounted arm pad, 330mm long, sliding. Standard option.

Armrest Height - adjustable from 18-26cm (2cm increments)

Armrest Pad Position - adjustable from 33-45cm (2cm increments)

Accessories

Pelvic Strap



Calf Strap



Headrest Kit



IV Pole



Meal Tray



Oxygen Bottle Carrier



Elevating Legrests



Left

Right

Both

Stump Support



Left

Right

Both

Requires Cushion

Yes Specify Cushion: _____

No

Notes