


# Staebel General Purpose Sling with Head Support


Product Code: 27GPH02

The Staebel General Purpose Head Support sling, as the name suggests, is a full body support style sling and is recommended for general lifting and in situations of lifting from the floor, and for the less secure patient.






**Mesh  
Available**



**Fabric  
Available**



**200 kg  
SWL**



### Specifications:

#### Mesh

Product Code	27GPHS02XS	27GPHS02S	27GPHS02M	27GPHS02L	27GPHS02XL	27GPHS02XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	820 mm	790 mm	920 mm	1110 mm	1270 mm	1360 mm
B	210 mm	230 mm	240 mm	260 mm	270 mm	290 mm
C	750 mm	780 mm	920 mm	1000 mm	1020 mm	1050 mm
D	1310 mm	1280 mm	1410 mm	1630 mm	1680 mm	1780 mm

#### Fabric

Product Code	27GPHS01XS	27GPHS01S	27GPHS01M	27GPHS01L	27GPHS01XL	27GPHS01XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	820 mm	790 mm	920 mm	1110 mm	1270 mm	1360 mm
B	210 mm	230 mm	240 mm	260 mm	270 mm	290 mm
C	750 mm	780 mm	920 mm	1000 mm	1020 mm	1050 mm
D	1310 mm	1280 mm	1410 mm	1630 mm	1680 mm	1780 mm