

# Staebel Slideezi Support Straps

Product Code: 210064

Support Straps used to reduce the risk of over-stretching and reaching in frequent supine transfers.





These Staebel Slideezi Support Straps are used to reduce the risk of over-stretching and reaching in frequent supine transfers.

A new concept in support straps, these can be wiped down and/or washed in soapy water, rinsed and then left to dry. Used to reduce the risk of over-stretching and reaching in frequent supine transfers, they simply clip onto the handles of the flat Slideeezi with handles and other patient handling equipment.

#### Product Features:

- For use in supine transfers
- Reducing the risk of overreaching and stretching
- Comes in a pair



#### Specifications:

Product Code	210064
Width	50 mm
Length	750 mm
SWL	100 kg