

# Staelbel Disposable Procedure Slide

Product Code: 210068

The Disposable Procedure Slide is a simple slide sheet that allows health professionals to easily move patients into the lithotomy position, ready for surgery.



## Product Features:

- Adjust patient position mid-surgery
- Slide user back up the operating table post-surgery
- Disposable slides help teams uphold excellent infection control
- Slippery material enables easy patient repositioning
- Integral locking pads ensure safety and security
- Can remain under the patient
- 25 per pack

## Specifications:

Product Code	210068
Width	410 mm
Length	410 mm