

Staelbel Disposable Rolling Slide

Product Code: 210134

The Disposable Procedure Slide is a simple slide sheet that allows health professionals to easily move patients into the lithotomy position, ready for surgery.



The Staebel Disposable Procedure Slide is a simple slide sheet that allows health professionals to easily move patients into the lithotomy position, ready for surgery. It is ideal for head-down procedures where the patient is positioned at an incline of up to 25 degrees and can be left under the patient.

Specifications:

Product Code	210134
Width	1000 mm
Length	2000 mm

Product Features:

- Slide a patient into a precise position (fine adjustments achievable), usually the Trendelenburg position
- Lock patient position with the sheet's unique and integral locking panels
- Position patients at inclines up to 25 degrees, facing head-down
- Slide user back up the operating table post-surgery
- Disposable slides help teams uphold excellent infection control
- Slippery material enables easy patient repositioning
- Integral locking pads ensure absolute safety and security
- 25 per pack