

# Stael Easy Dress Arms

Product Code: 210159

Designed to slide over a person's arm and ease the putting on of upper body garments.



A pair of long slipper sleeves that can slide over a person's arm and ease the putting on of upper body garments. For best practice: The sleeves work especially well with a flat slide sheet underneath the bed bound person to facilitate the slide.

## Product Features:

- Help carers put shirts/pyjamas on a bed bound or someone who requires a lot of assistance
- Easy wipe clean finish to promote infection control
- Slippery coating enables easy, dignified dressing

## Specifications:

Product Code	210159	210161
Size	Medium	Large
Width	280 mm	290 mm
Length	1070 mm	1090 mm