

Stael Anti Shear Mattress Cover

Product Code: 210290

A simple answer to the problem of constantly putting slide sheets in and out – lightweight and easy to use. The Anti Shear Mattress Cover's is simple as it makes it so much easier to adjust someone's position without disturbing them and reduces the risk of injury for staff.



Thinner than a bed sheet, this product fits over the mattress and does not impinge on any pressure reducing qualities the mattress may have. You can sit the client/patient up in bed without the risk of them slipping down, and it only becomes fully slippery when you introduce a flat slide sheet on top to perform the manoeuvre that needs to be performed.

Specifications:

Product Code	210290
Width	900 mm
Length	2000 mm

Product Features:

- Specialist technical fabric is fully breathable, lightweight and quiet
- Thin and light weight fabric
- Washable at 40°C to 74°C
- When using the base sheet to lie on, independent people can turn more easily in bed
- Breathable fabric means that people will not sweat or get too hot



Lightweight
Material



Easy to Use
Design