

# Staebel ComfortCare Bed Top Sheet

Product Code: 210295

The ComfortCare is an in-bed sliding system for the care of the long term or 24 hour patient.





### Product Features:

- Patients can remain on this indefinitely, either directly or under a cotton bed sheet
- Enhances patient care and clinical outcomes
- Polycotton system, with low-friction inner linings
- Two panel system enables up/down manoeuvres, transferring and turning into prone position
- Integral locking panel means safe position can be secured
- Suitable for single carers
- Washable at 40°C to 80°C and low temperature tumble dry
- Spare panels available

### Specifications:

**Product Code** 210295

**Width** 1400 mm

**Length** 2000 mm

**SWL** 200 kg

