

Staebel ComfortCare Bed Base Sheet

Product Code: 210296

The ComfortCare is an in-bed sliding system for the care of the long term or 24 hour patient.





Product Features:

- Allows carers to move a patient in all directions
- Patients can remain on the 4-way slide sheet indefinitely, either directly on top of the sheet or under a cotton bed sheet
- Carers can re-position patients without disturbing them or handling them directly
- Reduces risk of injury to carers
- 100% cotton system with low-friction inner linings
- Two-panel system enables up/down manoeuvres, turning, swiveling, sitting up or transferring
- Integral locking panels means safe position can be secured, even when sitting
- Washable at 40–74 °C and suitable for low-temperature tumble dry
- Spare panels available



Specifications:

Product Code	210296
Width	900 mm
Length	2000 mm
SWL	225 kg