


Staebel Bed Transfer Assist

Product Code: 210315

Transfer Assist Device designed to assist users in supporting the lower limbs during functional transfers.







The Staebel Bed Transfer Assist has been designed to assist in supporting the lower limbs during functional transfers. Simple and brilliant, our Bed Transfer Assistant allow one carer to safely raise a person's legs and swivel them around into bed and vice versa.

Product Features:

- Cotton upper for comfort
- Raising legs and swiveling in a sitting position by a single carer
- Ultra-slippery glide sheet makes transfers smooth and easy
- Padded cushion ensures patient comfort
- Rotary disc helps patients to rotate easily
- Washable at 74° C (glide sheet only)



Specifications:

Product Code	210315
Width	600 mm
Length	640 mm
Thickness	30 mm
SWL	100 kg