

Staebel Transfer Board Slide

Product Code: 210360

This device allows for a sitting transfer across the board without risk of shearing tissue damage.



The Transfer Board Slide is a low-cost device that allows users to independently transfer from a bed to a chair, wheelchair or car seat without risk of shearing-induced tissue damage. Simply lay the slide sheet on top of your transfer board, tuck it under your bottom and slide from one surface to another to enable easy transitioning between surfaces.

Simple and highly effective, this transfer board glide sheet is only as wide as the board itself and has a panel of anti-slip that grips the top of the board. This helps to prevent slippage during use, keeping patients safe and secure.

The slide sheet is suitable for constant, long-term use. Alternatively, the Patient-Specific Transfer Board Slide is often better-suited to patients whose reduced mobility is temporary, especially in high turnover areas.

Product Features:

- Can be folded up and popped in a bag
- Helping people move, both independently or with carer assistance
- Easily washable to promote excellent infection control

Specifications:

Product Code	210360
Width	300 mm
Length	700 mm
SWL	500 kg