

Staebel In Bed Turner Base Sheet

Product Code: 2103TBS

The InBed Turner Base Sheets allow safe and easy repositioning for people in bed enabling them to turn easily from side to side but not slip down.



The InBed Turner Base Sheet is part of our renowned satin range for in-bed repositioning and patient care. All InBed Turner Base Sheets allow safe and easy repositioning for people in bed enabling them to turn easily from side to side but not slip down.

Extra safety is provided by the striped cotton edges on all of our Base Sheets. They can be used on their own or combined with one of the Top Sheets for those who need more in bed repositioning care.



Product Features:

- Reduces risk of tissue and nerve damage
- Helps reduce pressure sores developing
- No training required
- Easy turning and repositioning with little effort
- Washable Design
- Full length or midi length satin
- Central satin panel is smooth for ease of movement
- Anti-slip cotton polyester on edges reduces risk of slipping
- Cotton polyester for enhanced comfort and repeat washing

Specifications:

Product Code	210321	210323	210325
Size	Single	King Single	Double
Width	900 mm	1050 mm	2300 mm
Length	2000 mm	2000 mm	2000 mm
SWL	300 kg	300 kg	300 kg