

Staebel Transfer Platform

Product Code: 212150

Designed as a platform for use by healthcare professionals to assist in patient transfers in a standing position over a short distance.





Staebel Transfer Platform is a manually-operated device designed as a platform for use by healthcare professionals to assist in the physical moving (transfer) of a patient in an upright (standing) position over a short distance. The device consists of a wheeled platform, a surrounding framework, kneepad and handlebar.

Product Features:

- Fast, stable and safe transport
- Assists clients with moving from bed to the chair
- Knee support, handlebar and safety belt to ensure safe transfers
- Safety Belt can be added on for additional support and security

Specifications:

Product Code	212150
Width	620 mm
Height	1180 mm
SWL	150 kg

