

# Oval Board Slide

Product Code: 2150B

Used to facilitate placing of leg supports during application of lifting slings with divided leg supports in a seated position.



Another popular product from the RoMedic range the EasyGlide black oval boards can be used to facilitate placing of leg supports during application of lifting slings with divided leg supports in a seated position. Due to their unique design and ability to flex these also eliminate the friction under the user and make it very easy to position the user further back in the chair.

## Product Features:

- Easy to Use Design
- Supplied in pairs

## Specifications:

Product Code	215040	215045
Width	850 mm	450 mm
Length	390 mm	200 mm
Thickness	3 mm	3 mm
SWL	135 kg	135 kg