


Foot Stool

Product Code: 218090

FootStool is used to raise the position of the user's knees to enable an easier and safer transfer in a seated position.



Product Features:

- For a sitting transfer from a wheelchair, FootStool is often used in combination with EasyGlide sliding board.
- An indentation in the front edge of FootStool enables the carer to stand as close as possible, to give better posture.
- In the rear edge there is a more shallow indentation for a toilet chair or other obstacles.
- Made of cellular plastic, a sturdy material which doesn't slide on the underlying surface.
- The surface is easy to clean with a soap or alcohol solution.

Specifications:

Product Code	218090
Width	380 mm
Length	320 mm
Thickness	100 mm