




ActiveCare General Purpose Sling with Head Support

Product Code: 227210

The General Purpose Head Support sling, as the name suggests, is a full body support style sling and is recommended for general lifting and in situations of lifting from the floor, and for the less secure patient. The loop style attachments have three lengths at each connection point and are used on a patient lifter with a yoke style attachment.




**Mesh
Available**



**Fabric
Available**



**Parasilk
Available**



Specifications:

Mesh

Product Code	227210	227126	227127	227128	227248	227212
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	795 mm	895 mm	995 mm	1120 mm	1240 mm	1360 mm
B	190 mm	210 mm	230 mm	245 mm	250 mm	250 mm
C	685 mm	790 mm	895 mm	995 mm	1080 mm	1125 mm
D	1085 mm	1270 mm	1440 mm	1595 mm	1735 mm	1860 mm
E	585 mm	650 mm	700 mm	740 mm	780 mm	810 mm

Fabric

Product Code	227209	227247	227067	227069	227205	227211
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	795 mm	895 mm	995 mm	1120 mm	1240 mm	1360 mm
B	190 mm	210 mm	230 mm	245 mm	250 mm	250 mm
C	685 mm	790 mm	895 mm	995 mm	1080 mm	1125 mm
D	1085 mm	1270 mm	1440 mm	1595 mm	1735 mm	1860 mm
E	585 mm	650 mm	700 mm	740 mm	780 mm	810 mm

Parasilk

Product Code	227410	227411	227412	227413	227414
Size	Extra Small	Small	Medium	Large	Extra Large
A	795 mm	895 mm	995 mm	1120 mm	1240 mm
B	190 mm	210 mm	230 mm	245 mm	250 mm
C	685 mm	790 mm	895 mm	995 mm	1080 mm
D	1085 mm	1270 mm	1440 mm	1595 mm	1735 mm
E	585 mm	650 mm	700 mm	740 mm	780 mm