

Digital Wheelchair Scale

Product Code: 2664

With its high capacity of 360 kg, this digital wheelchair scale makes it easy to weigh bariatric patients.



With its high capacity of 360 kg, this scale makes it easy to weigh bariatric patients. Thanks to its extremely low platform and the self-lowering ramps, this scale is easy to access. The integrated pre-TARE function determines in seconds the new weight of the patient utilizing internal weight memory. The swivelling LCD display with large letters is easy to read both for the patient and for the carer. The scale can be folded up in no time and can be transported anywhere thanks to its unique roller castors.

Product Features:

- Power pack is included
- Optional : RS232 adapter set
- 24 month warranty

Specifications:

Product Code	2664
Platform Width	760 mm
Platform Depth	930 mm
Platform Height	55 mm
Graduation	50 g
Device Weight	24.5 kg
SWL	360 kg