

Staebel Salsa 200 Standing Lifter

Product Code: 27700

Designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position.



Designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position. The Salsa 200 is ideal for transferring from bed to chair, or chair to commode. Allegro Concepts standing slings are manufactured to suit the Salsa, allowing superior access for dressing and personal hygiene.

With a small footprint and low overall weight, the Salsa 200 is highly maneuverable. A high capacity safe working load of 200kg means the Salsa is equally at home in an institutional or home environment.

Designed and tested to comply with:
AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

Product Features:

- Australian made quality product
- Electric leg adjustment
- Patient friendly lift arms enable multiple grip positions during lift
- Removable footplate for ease of transfers
- Height and angle adjustable contoured knee pad with removable calf strap
- Comes complete with battery pack and charger
- T.G.A. 136 065
- Safe Working Load of 200 kg

Specifications:

Product Code	27700
Base Length	1010 mm
External Base Width Closed	665 mm
External Base Width Open	1030 mm
Internal Base Width Closed	565 mm
Internal Base Width Open	930 mm
Base Width Adjustment	Electric
Under Bed Clearance	130 mm
Minimum Height Lift Arm	880 mm
Maximum Height Lift Arm	1685 mm
Under Base Height	55 mm
Lifter Weight	51 kg
Minimum Storage Height	1190 mm
Footplate	Removable
Kneepad	Height Adjustable
SWL	200 kg

