



# Staebel Pivot Frame Sling

Product Code: 27P02

The Pivot sling is designed for use with the pivot/tilt frame. It attaches with the keyhole tags engineered to the highest standards in strength and size tolerances and feature auto clip on action (no safety straps required). It also features sling centre marking and positioning handles to assist correct positioning and woven label that will not wash out or fade.



## Specifications:

### Mesh

Product Code	27P02XS	27P02S	27P02M	27P02L	27P02XL	27P02XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	720 mm	870 mm	1000 mm	1080 mm	1170 mm	1320 mm
B	190 mm	230 mm	270 mm	260 mm	270 mm	290 mm
C	580 mm	830 mm	920 mm	1010 mm	1090 mm	1180 mm
D	870 mm	1220 mm	1310 mm	1410 mm	1540 mm	1710 mm

### Fabric

Product Code	27P01XS	27P01S	27P01M	27P01L	27P01XL	27P01XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	720 mm	870 mm	1000 mm	1080 mm	1170 mm	1320 mm
B	190 mm	230 mm	270 mm	260 mm	270 mm	290 mm
C	580 mm	830 mm	920 mm	1010 mm	1090 mm	1180 mm
D	870 mm	1220 mm	1310 mm	1410 mm	1540 mm	1710 mm