


Stabel Pivot Care Frame Sling


Product Code: 27PC02

The Stabel Pivot Care Frame sling, is a comfortable, supportive sling with easily defined attachment points. Stabel Pivot Care Slings are used with hoists fitted with a pivot attachment. Sling provides easy removal of clothing for toilet transfers. Recommended for general toilet transfers for the less secure patients. Available in Mesh or Fabric and in sizes extra small to extra extra large.






Mesh Available



Fabric Available



200 kg SWL



Specifications:

Mesh

Product Code	27PC02XS	27PC02S	27PC02M	27PC02L	27PC02XL	27PC02XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	680 mm	840 mm	1020 mm	1030 mm	1220 mm	1240 mm
B	170 mm	210 mm	220 mm	220 mm	210 mm	210 mm
C	540 mm	630 mm	700 mm	790 mm	830 mm	790 mm
D	870 mm	1220 mm	1240 mm	1360 mm	1430 mm	1480 mm
T	980 mm	1180 mm	1220 mm	1320 mm	1440 mm	1660 mm

Fabric

Product Code	27PC01XS	27PC01S	27PC01L	27PC02L	27PC01XL	27PC01XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	680 mm	840 mm	1020 mm	1030 mm	1220 mm	1240 mm
B	170 mm	210 mm	220 mm	220 mm	210 mm	210 mm
C	540 mm	630 mm	700 mm	790 mm	830 mm	790 mm
D	870 mm	1220 mm	1240 mm	1360 mm	1430 mm	1480 mm
T	980 mm	1180 mm	1220 mm	1320 mm	1440 mm	1660 mm