


# Staebel Standing Loop Sling, Velcro Fastening

Product Code: 27STL04

The Staebel Standing Loop Sling is for use with the stand up hoists with loop attachment. It features optimum back support and under arm padding for a safe and comfortable lift and also allows for the easy removal of clothing for toilet transfers. It features a torso band for security and stability with Buckle or Velcro fastening.



## Specifications:

Product Code	27STL04S	27STL04M	27STL04L	27STL04XL	27STL04XXL
Size	Small	Medium	Large	Extra Large	Extra Extra Large
A	730 mm	880 mm	1060 mm	1120 mm	1430 mm
B	260 mm	290 mm	370 mm	410 mm	470 mm
C	900 mm	960 mm	1020 mm	1220 mm	1660 mm

Issue Date: 190825  
Revision Date: 190825 | V1