


NEOS 150 Patient Lifter

Product Code: 2880

A compact designed medium sized heavy duty mobile patient lifting hoist with a safe working load of 150kg ideally suited as a compact and lightweight hoist in a home care environment.



The NEOS 150 Patient Lifter is a compact ergonomic designed medium sized heavy duty mobile patient lifting hoist with a safe working load of 150kg. The NEOS 150 whilst being compact, provides superior patient space. The NEOS 150 in fact offers more patient space than most other larger size lifters. The swivel feature of the detachable yoke allows a patient to be rotated through 360° at full lift height, thus making the carer's job easier and safer.

Product Features:

- Folding Aluminium Lifter with a Safe Working Load (SWL) of 150kg
- Quick release folding structure enables lift to be conveniently folded for storage and transfer
- Electric raise and lower powered by Linak
- Foot pedal control for leg adjustment
- Equipped with 75mm front castors and 100mm locking rear castors
- Quick release yoke

Specifications:

Product Code	2878
Base Width Min	615 mm
Base Width Max	965 mm
Folded Overall Length	1230 mm
Folded Overall Width	400 mm
Lift Range - Yoke	410 - 1650 mm
Castor Size	75 / 100 mm
Under Bed Clearance	135 mm
Lifter Weight	35 kg
SWL	150 kg

Issue Date: 251022
Revision Date: 201025 | V2