

# Alto MKII 200 Patient Lifter

Product Code: 2891

The ultimate medium sized heavy duty mobile patient lifting hoist.



## The lifter of choice for institutional use or home care settings.

The Alto MKII 200 Patient Lifter is a medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design. The Alto MKII 200 Lifter, whilst being compact, provides superior patient lift space through the use of SPACELINK™ technology.

The Alto MKII 200 Lifter is the lifter of choice for institutional use or home care where a high capacity lift is required. Electric leg adjustment is a standard feature. It utilises the new SPACELINK™ concept, a breakthrough for optimised patient space, that enhances floor lift capability and increases max. lifting height. With maintenance free bushes in all the major moving parts this hoist is an Australian Made example of Quality and Strength that is hard to find in other hoists.

The Alto MKII 200 Patient Lifter is the lifter of choice for institutional use or home care where a high capacity lift is required. Designed and tested to comply with: AS ISO 10535-2002 Hoists for transfer of disabled persons - Requirements and test methods.

### Product Features:

- Electric leg adjustment is a standard feature
- Battery life indicator on lifter as well as handset
- Additional up and down buttons on control box
- Emergency stop button
- Supplied complete with battery, battery charger and instruction manual
- Australian Owned and Made
- Compatible with our Staebel Slings & a range of hoist attachments



Standard Yoke



4 Point Yoke



Tri Hook Yoke



Pivot Frame



Powered Pivot Frame



Integrated Weigh Scale

### Specifications:

<b>Product Code</b>	2891
<b>Base Length</b>	1150 mm
<b>External Base Width Closed</b>	700 mm
<b>External Base Width Open</b>	1130 mm
<b>Internal Base Width Closed</b>	600 mm
<b>Internal Base Width Open</b>	1030 mm
<b>Base Width Adjustment</b>	Electric
<b>Underbed Clearance</b>	130 mm
<b>Minimum Height Yoke</b>	300 mm

<b>Maximum Height Yoke</b>	1920 mm
<b>Average Lift Range</b>	1620 mm
<b>Minimum Storage Height</b>	1350 mm
<b>Under Bed Height</b>	55 mm
<b>Hoist Reach</b>	610 - 1060 mm
<b>Hoist Reach at Maximum Height</b>	515 mm
<b>Mass of Heaviest Part</b> (Base Assembly)	24 kg
<b>Lifter Weight</b>	46 kg
<b>SWL</b>	200 kg

Issue Date: 280415  
Revision Date: 201025 | V4