

Staebel Alto MKII 200 Patient Lifter

Product Code: 2891

The ultimate medium sized heavy duty mobile patient lifting hoist.



The lifter of choice for institutional use or home care settings.

The Staebel Alto MKII 200 Patient Lifter is a medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design. The Alto MKII 200 Lifter, whilst being compact, provides superior patient lift space through the use of SPACELINK™ technology.

The Alto MKII 200 Lifter is the lifter of choice for institutional use or home care where a high capacity lift is required. Electric leg adjustment is a standard feature. It utilises the new SPACELINK™ concept, a breakthrough for optimised patient space, that enhances floor lift capability and increases max. lifting height. With maintenance free bushes in all the major moving parts this hoist is an Australian Made example of Quality and Strength that is hard to find in other hoists.

The Staebel Alto MKII 200 Patient Lifter is the lifter of choice for institutional use or home care where a high capacity lift is required. Designed and tested to comply with: AS ISO 10535-2002 Hoists for transfer of disabled persons - Requirements and test methods.

Product Features:

- Electric leg adjustment is a standard feature
- Battery life indicator on lifter as well as handset
- Additional up and down buttons on control box
- Emergency stop button
- Supplied complete with battery, battery charger and instruction manual
- Australian Owned and Made
- Compatible with our Staebel Slings & a range of hoist attachments



Standard Yoke



4 Point Yoke



Tri Hook Yoke



Pivot Frame



Powered Pivot Frame



Integrated Weigh Scale

Specifications:

Product Code	2891	Maximum Height Yoke	1920 mm
Base Length	1150 mm	Average Lift Range	1620 mm
External Base Width Closed	700 mm	Minimum Storage Height	1350 mm
External Base Width Open	1130 mm	Under Bed Height	55 mm
Internal Base Width Closed	600 mm	Hoist Reach	610 - 1060 mm
Internal Base Width Open	1030 mm	Hoist Reach at Maximum Height	515 mm
Base Width Adjustment	Electric	Mass of Heaviest Part (Base Assembly)	24 kg
Underbed Clearance	130 mm	Lifter Weight	46 kg
Minimum Height Yoke	300 mm	SWL	200 kg