


Transfer Board

Product Code: 292TB

A semi-rigid board constructed from a lightweight, durable, smooth plastic this product is used to transfer a patient in the lying down position from surface to surface.



Product Features:

- Is easily cleaned by wiping over with mild detergent
- Radio-lucient
- Rounded edges and radius corners facilitate a smooth transfer
- Underside features grip sections
- Can be mounted on a hook on wall
- 150kg SWL

Specifications:

Product Code	2923	2924
Width	150 mm	180 mm
Length	580 mm	580 mm
SWL	135 kg	135 kg