

Seca Digital Chair Scale

Product Code: 2954

Consideration and comfort – these are particularly prominent aspects of this extremely stable model with its circular steel tube frame.



A comfortable seat provides a safe position, the swivelling armrests make it easy to move the patient and the adjustable footrests increase ease of operation. The locking rear wheels ensure that it is safe to sit down and stand up at all times. The position of the LCD display unit behind the chair makes it easy to read and has an In-built recalibration function.

Product Features:

- Flip up armrests and footrests
- 24 month warranty
- Functions : Automatic switch off, TARE, HOLD, BMI, kg/lbs/ sts switch over

Specifications:

Product Code	2954
Width	563 mm
Depth	978 mm
Graduation (>150kg)	50 g
Graduation (<150kg)	100 g
Overall Height	330 mm
SWL	300 kg