

Staebel Chair Scale

Product Code: 2957

The Staebel Chair Scale with its patented design maintains reliability and uncompromising precision.



The Staebel Chair Scale with its patented design maintains reliability and uncompromising precision. The compact design ensures the most cost efficient shipping method and reduces the risk of carer strain particularly when weighing an obese/bariatric patient. The unique base frame design gives care staff assurance that scale is compatible with other equipment utilised in care setting.

Product Benefits:

- Twin wheel castors are standard and do not create excessive noise when wheeling scale along a corridor
- One piece footplate flips up and has been designed with skin tear prone patients in mind
- Armrests flip back for scale to be used in a bariatric care setting
- Base frame comes standard with a level indicator
- BMI function/ Tare/Zero/Hold/Print/kg/lb
- Scroll Down/Up on LCD screen
- 5 digit ¾” LCD screen that can be tilted to accommodate different care staff heights
- Automatic re-weigh
- Up-datable firmware
- Ports: RS232 and USB
- Lightweight yet durable construction - only weighs 17.5 kg
- 24 month warranty

Specifications:

Product Code	2957
Overall Width	630 mm
Overall Height	965 mm
Seat Width	430 mm
Seat Depth	420 mm
Weight	17.5 kg
SWL	300 kg

