

C Series C450 Power Traverse Ceiling Hoist

Product Code: 7323149

The C450 has been designed with a full set of standard features and can be safely operated by a single caregiver. Built for quiet operation, the hoist provides a smooth & effortless transfer.



Quiet and smooth, our compact C-Series ceiling lifts feature a pneumatic hand control that has no electrical components for additional safety. The large buttons make it easy to operate the lift with clear directional indicators.

Product Features:

- Smallest Hoist of its class which allows for 180mm of additional lifting height
- Quick lifting speed
- Emergency pull cord
- Modular design for simplified maintenance
- Quick charge batteries for reduced service time
- Audible low battery indicator & visual battery/ charge level display
- Digital readout
- Auto shut off when not in use, to save battery life
- 2.2 strap vertical movement allows you to easily lift from the floor
- Pneumatic Hand Controls – 100% waterproof
- Manual option – for emergency raising or lowering
- Available with return to charge option

Specifications:

Product Code	7323149
Total Length	318 mm
Total Width	218 mm
Lifting Range	2082 mm
Weight	10 kg
SWL	204 kg



Easy to Use Design



204kg SWL