

Safe Sense 2 Stand-Alone Bed Exit System

Product Code: 50-0984

Sustainable fall prevention for utmost safety. This is what the SafeSense® bed exit assistance system of wissner-bosserhoff stands for.



SafeSense® signals the nursing staff that the resident has left the bed or is about to leave it. This allows falls to be detected earlier than through routine checks. Essentially, the bed-exit assistance system works like a traffic light.

SafeSense® can considerably increase the level of safety, especially during night and weekend shifts with fewer staff. Especially at night, the bed exit assistant system can minimize the need for side guards. The freedom of the resident is upheld. An automatic night light helps the resident to safely find the way to the bathroom at night. It is automatically switched on when the resident leaves the bed and is switched off on return. The interval of the nurse call timer can be set flexibly according to the individual degrees of mobility of the residents.

Product Features:

Greater safety – fewer falls/restraints

- Faster response time in the event of falls
- Greater safety during night and weekend shifts
- Ideal in combination with sentida and image 3 low nursing beds and SafeFree® side guards
- Systematic implementation of the findings from the initiative for the reduction of restraints, the "Werdenfelser Weg" (Werdenfels care pathway)

Advantages for your Aged Care Home or Health Care Facility

- Investment security
- Resident safety in any working position
- No additional costs for safety equipment
- Efficient use of restraining belts

Individual Protection

- Sensor-switchable night light enables light e.g. on the way to the toilet at night
- The nurse call timing is flexibly set from 0 to 30 min. according to the mobility of the resident

Compatibility for retrofitting

- On all wissner-bosserhoff beds
- On various third-party bed

Stand Alone Version

- Available as a cable or wireless solution.
- This product can be used as an added option or retrofit kit for wissner-bosserhoff beds movita sc and carisma sc, as well as for various brand of beds

A trial maybe necessary to ensure compatibility



Specifications:

Product Code 50-0984

Length 800 mm

Also require adapter cable