

Cura1 Premium Monitor

Product Code: 92611

The workhorse of your facility.
Compatible with Premium, SafeTime and
Universal (if corded) pads.



Product Features:

- Change Pad (SafeTime) indication*
- Pause mode to attend to a patient without triggering the alarm
- Connects with nurse call system
- Tamper-proof reset option using magnetic caregiver key
- Choice of tone, and volume control for alarm
- 0, 2 sec, and 15 min. alarm delay settings
- Port to connect with nurse call system (NC/NO).
- Operates on batteries and / or optional external plug pack (2406)
- Low battery notification
- Input Lost notification (pad lost, nurse call, AC lost)
- Comes with mounting bracket, hanger strap, soft cover, reset key and batteries
- Warranty: 24 months

* SafeTime feature only applicable when used with SafeTime bed and chair pads.

Specifications:

Product Code	92611
Length	120 mm
Width	77 mm
Height	34 mm
Operating Temperature	5°C - 40°C
Battery Life	> 1 year
Weight (excluding batteries)	140 g