

Cura1 Pull String Bed Monitor with Magnet

Product Code: 92617

Simple but very effective pull string activates alarm. Simply clip the pull-string to patient's clothing to raise an alarm when they get out of their chair or bed.



Product Features:

- Pull string activates alarm
- Affordable price
- Alarm memory and auto re-arm
- Battery operated (battery included)
- Adjustable alarm volume
- Red LED alarm indicator light
- Low battery detection and indication
- Warranty: 12 months

Specifications:

Product Code	92617
Length	120 mm25
Width	77 mm
Height	120 mm