

Cura1 Cordless Chair Pad with Transmitter

Product Code: 92661

Cura1 Cordless Chair Pad with Transmitter to raise an alarm when someone gets off their chair.



Product Features:

- Simple operation and maintenance
- Includes transmitter
- For use with Cura1 cordless monitor
- Easy to clean
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- 12 month life expectancy

Must be used with 92700 Monitor

Specifications:

Product Code	92661
Length	255 mm
Width	320 mm
SWL	30 - 200 kg