

Cura1 Cordless Bed Pad with Transmitter

Product Code: 92662

Cura1 Cordless Bed Pad with Transmitter to raise an alarm when someone gets out of their bed.



Product Features:

- Pad includes pouch that contains a transmitter
- Guaranteed 1 year life expectancy
- Easy to clean
- Non-allergenic soft PVC
- For use with cura1 cordless bed monitors
- 50 meter range (*Subject to environmental conditions)
- Vent holes in the neck reduce the risk of moisture

Must be used with 92700 Monitor

Specifications:

Product Code	92662
Length	510 mm
Width	760 mm
SWL	30 - 200 kg