

# Cura1 Cordless Floor Pad with Transmitter

Product Code: 92663K

To be positioned beside a users bed to alert staff when they put their feet on the floor, attempting to leave the bed.



## Product Features:

- Simple operation
- Easy to clean
- Virtually no false alarms over the lifetime of the pad
- Made from heavy duty non-slip material
- Generous size compared to others on the market
- Quality Australian made product
- 12 month life expectancy
- 50 meter range (\*Subject to environmental conditions)
- Includes transmitter

*Must be used with 92700 Monitor*

## Specifications:

Product Code	92663K
Length	1220 mm
Width	610 mm
SWL	30 - 200 kg