

# Cura1 Six Channel Cordless Monitor

Product Code: 92700

No more cords between monitor and pads. Can be used with up to six Cordless pads, mats or other triggering devices.



## Product Features:

- Modern ergonomic design
- Monitor up to 7 devices simultaneously - a single Premium pad / mat and up to 6 Cordless devices.
- The 2700 comes with batteries, durable ABS plastic enclosure and Caregiver Reset Key
- Alarm memory and auto re-arm
- Amber LED visual alarm indicator
- Caregiver secure reset option
- Battery operated or optional external plug pack
- Up to 6 months battery life (nominal)
- Easy to access battery compartment
- Can be connected to a nurse call system
- 3 selectable alarm tones
- Easy thumbwheel adjustable alarm volume
- "Pause" feature allows a caregiver to attend to the patient without setting the alarm off
- Alarm Delay options
- Fault supervision for pressure pads (low battery, device lost)
- Alerts for lost nurse call connection, lost AC power supply & low battery
- Warranty: 24 months

## Safety Features:

- Non resettable "Pad Signal Lost" warning
- Caregiver secure reset option
- All safety alerts sound even if alarm volume set to silent, including low battery, AC lost, pad signal lost, change pad (SafeTime) & nurse call lost
- Cable strain relief & breakaway cord technology prevents damage & helps reduce risk of entanglement when used with Premium Pads & Mats

## Specifications:

Product Code	92700
Length	40 mm
Width	77 mm
Height	120 mm
Weight (without batteries)	150 g