



iBeam Patient Monitor

Product Code: 927413

Its narrow detection field and portability makes the 2721 iBeam one of our most versatile and popular devices in busy care facilities or at home.



Product Features:

- Simple operation and maintenance at an affordable price
- Narrow detection zone means a more precise target can be set. There are less false alarms in busy areas compared to conventional detectors
- Can be installed adjacent to a bed, in front of a chair or across a doorway, stairs, oven, heater or any potentially hazardous area
- Mounting bracket, screws and double-sided tape included
- On / Off switch to conserve battery life when not in use
- Low battery indicator
- Works with other cura1 cordless receivers or as a stand-alone alarm
- Accepts plug pack (2404, optional)
- Dedicated stand available (2723) for increased usability
- Warranty: 24 months

Specifications:

Product Code	92719
Width	78 mm
Depth	38.5 mm
Height	81.5 mm
Wireless Range	30 - 50 m*
Battery Life	6 months
Weight	100 g (excluding batteries)

*Actual transmission and detection distances may vary and be dependent on specific environmental conditions.