

# Safe Sense 2 Bed Exit System for Wi-Bo Beds

Designed to be a sustainable falls prevention system for utmost safety & protection.





SafeSense® signals the nursing staff that the resident has left the bed or is about to leave it. This allows falls to be detected earlier than through routine checks. Essentially, the bed-exit assistance system works like a traffic light.

SafeSense® can considerably increase the level of safety, especially during night and weekend shifts with fewer staff. Especially at night, the bed exit assistant system can minimize the need for side guards. The freedom of the resident is upheld. An automatic night light helps the resident to safely find the way to the bathroom at night. It is automatically switched on when the resident leaves the bed and is switched off on return. The interval of the nurse call timer can be set flexibly according to the individual degrees of mobility of the residents.

### Product Features:

#### Greater safety –fewer falls/restraints:

- Faster response time in the event of falls
- Greater safety during night and weekend shifts
- Ideal in combination with sentida and image 3 low nursing beds and SafeFree® side guards

#### Advantages for your Aged Care Home or Health Care Facility:

- Investment security
- Resident safety in any working position
- Legal certainty without having to comply with additional requirements
- No additional costs for safety equipment

#### Individual Protection:

- Sensor-switchable night light enables light e.g. on the way to the toilet at night
- The nurse call timing is flexibly set from 0 to 30 min. according to the mobility of the resident

#### Compatibility for retrofitting:

- On all wissner-bosserhoff beds
- On various brand of beds

#### Wi-Bo Bed Compatible:

- Integrated as an option
- As retrofit kit for all electrical nursing home beds by wissner-bosserhoff (manufactured in 2002 or later, exception: movita sc & carisma sc)
- Available with cable or as a wireless version

### Specifications:

Product Code	77025	77026	77015
Type	Wired for Sentida SC	Wireless for Sentida SC	Wireless for Duo9 Motors (> 2023 Beds)
Length	800 mm	800 mm	800 mm

