

CuroCell Static Air Cushion (S.A.C)

Product Code: 341SAC

The CuroCell® Static Air Cushion (S.A.C) is a static air cushion designed for comfort and pressure injury prevention.



160kg
SWL



Very High Risk
Cushion

2

Suitable for
Stage 2



Folds with the
Wheelchair



Ultra Lightweight
Design



Microclimate
Management

The CuroCell® Static Air Cushion (S.A.C.) is a static air cushion designed for comfort, pressure injury prevention, and as a support aid in the treatment of stage 1-2 pressure injuries.

Pre-inflated for immediate use, it requires no manual setup, making it simple and efficient for aged care, acute and home care environments. Constructed from a soft, durable material, it offers an optimal balance of comfort, support, and long-term performance.

Product Features:

- Air-venting perforations promote airflow and moisture dispersion to help maintain a healthy microclimate
- Soft, durable material supports long-term use without compromising comfort
- Lightweight, low-profile design is easy to position, handle, and clean
- Pre-inflated and patient-ready, eliminating the need for manual setup
- Reduces interface pressure to assist in the prevention and management of early-stage pressure injuries
- Low profile
- WARRANTY: 1 Year



Specifications:

Product Code	3417	3418
Width	450 mm (18")	500 mm (20")
Length	450 mm (18")	500 mm (20")
Thickness	50 mm (1.9")	50 mm (1.9")
SWL	160 kg	160 kg