

Staebel Relax Duo Profile Cushion

Product Code: 353DPC

Designed to quickly immerse into the cushion, while conforming to the users shape.




Designed to quickly immerse into the cushion, while conforming to the users shape.

The Staebel Relax Duo-Profile Cushion offers a pre-contoured surface giving excellent hip and leg alignment. The three layered design creates comfort for users with up to a medium pressure injury risk. The Duo-Profile uses a collapsible gel topper. This unique design allows you to quickly immerse into the cushion, while it conforms to your shape.

Product Features:

- Pre-contoured cushion assists in hip and leg alignment.
- Allows the user to immerse and envelop into the cushion, the gel topper conforms to your shape.
- Clear protective (PU) film coating with no seams is easy to wipe clean and is suitable for incontinent users
- Slow foam base and a memory foam mid layer.
- Suitable for those who are high risk of pressure injuries
- Safe Working Load of 130kg
- Includes a black waterproof wipe down outer cover
- 3 Layered Design:
 - Base is Fast Recovery Foam
 - Mid-section is Slow Foam/ Memory Foam
 - Top layer is medical-grade silicone



Specifications:

Product Code	35380	35381	35382	35383
Width	400 mm	430 mm	450 mm	450 mm
Length	400 mm	430 mm	400 mm	450 mm
Thickness	80 mm	80 mm	80 mm	80 mm
SWL	130 kg	130 kg	130 kg	130 kg