


# Staebel Zoncecare Concave Mattress

Product Code: 3-CCS

The Staebel Zoncecare Concave Medical Mattress features concave sides, encouraging the patient to lie centrally and thereby minimising the possibility of falling. It facilitates pressure distribution and has a vapour-permeable cover. Additionally, the heavy-duty base reduces the chance of wear due to sliding. It has designed to deliver a cost-effective solution across healthcare environments.



### Product Features:

- Deluxe Grade foam core
- Deluxe Grade Concave Sides encourage patients to lie in the centre of the mattress and help prevent fall risks
- 2-way stretch cover assists with pressure redistribution
- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Heavy-duty PVC base provides a durable base to promote cover longevity and minimise wear caused by the mattress sliding on the bed base
- Waterfall zipper cover to minimise fluid penetration and infection control risks

### Specifications:

Product Code	3-CCS-S	3-CCS-KS
Size	Single	King Single
Width	860 mm	1050 mm
Length	1960 mm	1960 mm
Thickness	150 mm	150 mm
SWL	150 kg	150 kg



150kg  
SWL



Low Risk  
Mattress



Prevention  
Mattress



Reduces  
Falls Risks