


Staebel Zoncecare Concave Mattress

Product Code: 3-CCS

The Staebel Zoncecare Concave Medical Mattress features concave sides, encouraging the patient to lie centrally and thereby minimising the possibility of falling. It facilitates pressure distribution and has a vapour-permeable cover. Additionally, the heavy-duty base reduces the chance of wear due to sliding. It has designed to deliver a cost-effective solution across healthcare environments.



Product Features:

- Deluxe Grade foam core
- Deluxe Grade Concave Sides encourage patients to lie in the centre of the mattress and help prevent fall risks
- 2-way stretch cover assists with pressure redistribution
- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Heavy-duty PVC base provides a durable base to promote cover longevity and minimise wear caused by the mattress sliding on the bed base
- Waterfall zipper cover to minimise fluid penetration and infection control risks
- WARRANTY: Foam core – 5 years , Cover – 1 years

Specifications:

Product Code	3-CCS-S	3-CCS-KS
Size	Single	King Single
Width	860 mm	1050 mm
Length	1960 mm	1960 mm
Thickness	150 mm	150 mm
SWL	150 kg	150 kg

Issue Date: 100325
Revision Date: 200126 | V3



150kg
SWL



Low Risk
Mattress



Prevention
Mattress



Reduces
Falls Risks