

Staebel Total Flow Plus Bariatric Active Air Mattress

Product Code: 3202301

Staebel Total Flow Bariatric Plus Mattress is a full replacement alternating mattress system designed as a cost-efficient solution for effective prevention of bariatric users at very high risk of developing pressure injuries.



Product Features:

- Nylon TPU independently removable cells
- Cell-on-cell design to prevent bottoming out and increase stability when the alternating mode is active
- Remote control to enable independence
- Pillow function/static head cells
- CPR quick release tag & anchoring straps
- Multi-stretch, antimicrobial, machine washable, waterproof & vapour permeable PU cover
- Micro-air loss cells in the torso region to support microclimate management
- Transport Mode
- Warranty: Mattress - 12 months, Pump - 24 months
- Standard: ISO 9001:2015 | ISO 13485:2016

Pump Features:

- Cycle Time: 10, 15, 20, 25 minutes
- Control panel lock to prevent unauthorised or accidental tampering
- 20-minute automatic timer for max inflate mode to reduce risk of prolonged max inflation
- Memory function in the case of a power outage
- Air Output: 8 litres/minute
- Power Input: AC220-240V, 50-60Hz
- Power Consumption: 7W
- Fitted with hanging brackets
- Dimensions (WxLxH): 240 x 345 x 125 mm



350kg
SWL



Very High Risk
Mattress



Prevention
Mattress



Bariatric
Use



Micro Air
Loss

A hand is holding a black cloth with two black straps attached, positioned below a row of blue and white ovals containing the letters 'c', 'b', and 'A'. The ovals are arranged in a row, with the letters 'c', 'b', and 'A' appearing in a sequence. The background is white.

Control Panel:

